


Product Datasheet

ORF37 (VZV) antibody (orb1101143)

Description	Mouse monoclonal antibody to ORF37 (VZV)
Species/Host	Mouse
Reactivity	Virus
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	The immunogen was VZV ORF37 protein.
Form/Appearance	Liquid
Concentration	1 mg/mL
Storage	Long term -20 °C, short term +4 °C. Avoid freeze-thaw cycles.
Note	For research use only
Isotype	IgG2b
Clonality	Monoclonal
Purity	Affinity purified
Expiration Date	12 months from date of receipt.



WB analysis of cells:
a=uninfected cells...