


## Product Datasheet

### ORF57 (VZV) antibody (orb1101154)

<b>Description</b>	Mouse monoclonal antibody to ORF57 (VZV)
<b>Species/Host</b>	Mouse
<b>Reactivity</b>	Virus
<b>Conjugation</b>	Unconjugated
<b>Tested Applications</b>	ELISA, WB
<b>Immunogen</b>	The immunogen consisted of full length ORF57 and was produced in E.coli
<b>Form/Appearance</b>	Liquid
<b>Concentration</b>	1 mg/mL
<b>Storage</b>	Long term -20 °C, short term +4 °C. Avoid freeze-thaw cycles.
<b>Note</b>	For research use only
<b>Isotype</b>	IgG2b
<b>Clonality</b>	Monoclonal
<b>Purity</b>	Affinity purified
<b>Expiration Date</b>	12 months from date of receipt.



WB analysis of cells:  
a=uninfected cells...