


Product Datasheet

ORF67 (VZV) antibody (orb1101158)

Description	Mouse monoclonal antibody to ORF67 (VZV)
Species/Host	Mouse
Reactivity	Virus
Conjugation	Unconjugated
Tested Applications	ELISA, IF, WB
Immunogen	E.coli produced ORF67 C-terminal domain (amino acids 294-354)
Form/Appearance	Liquid
Concentration	1 mg/mL
Storage	Long term -20 °C, short term +4 °C. Avoid freeze-thaw cycles.
Note	For research use only
Isotype	IgG2b
Clonality	Monoclonal
Purity	Affinity purified
Expiration Date	12 months from date of receipt.



WB analysis of cells:
a=uninfected cells...