


## Product Datasheet

### M57 (MCMV) antibody (orb1101168)

|                            |                                                               |
|----------------------------|---------------------------------------------------------------|
| <b>Description</b>         | Mouse monoclonal antibody to M57 (MCMV)                       |
| <b>Species/Host</b>        | Mouse                                                         |
| <b>Reactivity</b>          | Virus                                                         |
| <b>Conjugation</b>         | Unconjugated                                                  |
| <b>Tested Applications</b> | FACS, WB                                                      |
| <b>Immunogen</b>           | MCMV virions                                                  |
| <b>Preservatives</b>       | Liquid                                                        |
| <b>Concentration</b>       | 1 mg/mL                                                       |
| <b>Storage</b>             | Long term -20 °C, short term +4 °C. Avoid freeze-thaw cycles. |
| <b>Note</b>                | For research use only                                         |
| <b>Isotype</b>             | IgG1 kappa                                                    |
| <b>Clonality</b>           | Monoclonal                                                    |
| <b>Purity</b>              | Affinity purified                                             |
| <b>MW</b>                  | 48 kDa                                                        |
| <b>Expiration Date</b>     | 12 months from date of receipt.                               |



Immunofluorescence  
analysis of MEF-K  
cel...