


Product Datasheet

UL89 (HCMV) antibody (orb1101170)

Description	Mouse monoclonal antibody to UL89 (HCMV)
Species/Host	Mouse
Reactivity	Virus
Conjugation	Unconjugated
Tested Applications	ELISA, IF, WB
Immunogen	The immunogen consisted of UL89 N terminal domain and was produced in E.coli
Form/Appearance	Liquid
Concentration	1 mg/mL
Storage	Long term -20 °C, short term +4 °C. Avoid freeze-thaw cycles.
Note	For research use only
Isotype	IgG2a
Clonality	Monoclonal
Purity	Affinity purified
Expiration Date	12 months from date of receipt.



UL89
(HCMV)
antibody