
Product Datasheet

Human Angiopoietin-2 / ANGPT2 Protein (orb1101173)

Description

Human Angiopoietin-2, His Tag is

(HEK293). It contains AA Tyr 19 - Ph...

Reactivity

Human

Endotoxins

1.0 EU per µg

Conjugation

Unconjugated

Target

Angiopoietin-2

Preservatives

20 mM MOPS, 150 mM NaCl, 5 mM CHAPS, pH7.5

Form/Appearance

Powder

Storage

-20°C

Tag

C-10×His

Note

For research use only

Application notes

This protein carries a polyhistidine tag at the C-terminus. The protein has a calculated MW of 56.8 kDa. The protein migrates as 65-80 kDa under reducing (R) condition (SDS-PAGE) due to glycosylation.

Protein Sequence

AAI26201.1

Purity

95%

Source

Human Angiopoietin-2, His Tag (orb1101173) is expressed from human 293 cells (HEK293). It contains AA Tyr 19 - Phe 496 (Accession # O15123-1).

MW

56.8 kDa

Uniprot ID
[O15123](#)
NCBI
[AAI26201.1](#)
Expiration Date

6 months from date of receipt.