
Product Datasheet

Mouse Complement C3 Protein (orb1101180)

Description

Mouse Complement C3, His Tag is

(HEK293). It contains AA Ile 25 - Asn...

Reactivity	Mouse
Endotoxins	1.0 EU per µg
Conjugation	Unconjugated
Target	Complement C3
Preservatives	PBS, pH7.4
Form/Appearance	Powder
Storage	-20°C to -70°C for 12 months in lyophilized state; -70°C for 3 months under sterile conditions after reconstitution. For long term storage, the product should be stored at lyophilized state at -20°C or lower.
Tag	C-10×His
Note	For research use only
Application notes	This protein carries a polyhistidine tag at the C-terminus. The protein has a calculated MW of 186.1 kDa. The protein migrates as 70 kDa, 100-120 kDa and 160-180 kDa under reducing (R) condition (SDS-PAGE) due to glycosylation.
Purity	90%
Source	Mouse Complement C3, His Tag (orb1101180) is expressed from human 293 cells (HEK293). It contains AA Ile 25 - Asn 1663 (Accession # P01027-1).
MW	186.1 kDa
Uniprot ID	P01027
NCBI	NP_033908.2
Expiration Date	6 months from date of receipt.