

## Product Datasheet

# IRAK1 (phospho-Thr387) antibody (HRP) (orb110265)

|                            |  |
|----------------------------|--|
| <b>Description</b>         | IRAK1 (phospho-Thr387) antibody (HRP)  |
| <b>Species/Host</b>        | Rabbit   |
| <b>Reactivity</b>          | Human, Mouse, Rat  |
| <b>Conjugation</b>         | HRP  |
| <b>Tested Applications</b> | ELISA, IHC-Fr, IHC-P, WB   |
| <b>Immunogen</b>           | KLH conjugated Synthesised phosphopeptide derived from human IRAK1 around the phosphorylation site of Thr387 |
| <b>Target</b>              | Phospho-IRAK1 (Thr387)   |
| <b>Preservatives</b>       | 0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.  |
| <b>Form/Appearance</b>     | Liquid   |
| <b>Concentration</b>       | 1mg/ml   |
| <b>Storage</b>             | Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.                             |
| <b>Note</b>                | For research use only  |
| <b>Isotype</b>             | IgG  |
| <b>Clonality</b>           | Polyclonal   |
| <b>Uniprot ID</b>          | <a href="#">Q62406</a>   |
| <b>Dilution Range</b>      | WB=1:500-2000,IHC-P=1:100-500,IHC-F=1:100-500,ELISA=1:5000-10000   |
| <b>Expiration Date</b>     | 12 months from date of receipt.  |