

## Product Datasheet

# ROCK1 (phospho-Thr455/Ser456) antibody (HRP) (orb110401)

<b>Description</b>	ROCK1 (phospho-Thr455/Ser456) antibody (HRP)
<b>Species/Host</b>	Rabbit
<b>Reactivity</b>	Human, Mouse
<b>Conjugation</b>	HRP
<b>Tested Applications</b>	ELISA, IHC-Fr, IHC-P, WB
<b>Immunogen</b>	KLH conjugated synthesised phosphopeptide derived from human ROCK1 around the phosphorylation site of Thr455+Ser456
<b>Target</b>	phospho-ROCK1 (Thr455+Ser456)
<b>Preservatives</b>	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
<b>Form/Appearance</b>	Liquid
<b>Concentration</b>	1mg/ml
<b>Storage</b>	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
<b>Note</b>	For research use only
<b>Isotype</b>	IgG
<b>Clonality</b>	Polyclonal
<b>Uniprot ID</b>	<a href="#">Q13464</a>
<b>Dilution Range</b>	WB=1:500-2000,IHC-P=1:100-500,IHC-F=1:100-500,ELISA=1:5000-10000
<b>Expiration Date</b>	12 months from date of receipt.