

## Product Datasheet

### NOX2 antibody (orb1105610)

|                            |   |
|----------------------------|---|
| <b>Description</b>         | Rabbit polyclonal antibody to NOX2  |
| <b>Species/Host</b>        | Rabbit  |
| <b>Reactivity</b>          | Human, Mouse, Rat   |
| <b>Conjugation</b>         | Unconjugated  |
| <b>Tested Applications</b> | ELISA, ICC, IHC, WB   |
| <b>Target</b>              | NOX2  |
| <b>Preservatives</b>       | PBS   |
| <b>Form/Appearance</b>     | Liquid  |
| <b>Concentration</b>       | 1mg/ml  |
| <b>Storage</b>             | Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles. |
| <b>Note</b>                | For research use only   |
| <b>Clonality</b>           | Polyclonal  |
| <b>Expiration Date</b>     | 12 months from date of receipt.   |