

Product Datasheet

MRP4/Abcc4 antibody (orb11077)

Description	MRP4/Abcc4 antibody
Species/Host	Rabbit
Conjugation	Unconjugated
Tested Applications	ELISA, IHC-P
Immunogen	KLH conjugated synthetic peptide derived from human MRP4 (751-880/1325 aa)
Target	MRP4/ABCC4
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	149 kDa
Uniprot ID	O15439
Expiration Date	12 months from date of receipt.