

## Product Datasheet

# Anti-HEL Mouse IgG2b isotype control (orb1108554)

<b>Description</b>	Anti-HEL Mouse IgG2b isotype control
<b>Species/Host</b>	Mouse
<b>Conjugation</b>	Unconjugated
<b>Tested Applications</b>	Control
<b>Preservatives</b>	PBS pH 7.4
<b>Storage</b>	Store at -20...-70°C upon receipt. Avoid multiple freeze-thaw cycles
<b>Note</b>	For research use only
<b>Isotype</b>	mIgG2b
<b>Clonality</b>	Recombinant
<b>Uniprot ID</b>	NA
<b>Expiration Date</b>	12 months from date of receipt.