

Product Datasheet

Anti-HEL Mouse IgG1 isotype control (orb1108588)

Description	Anti-HEL Mouse IgG1 isotype control
Species/Host	Mouse
Conjugation	Unconjugated
Tested Applications	Control
Preservatives	PBS pH 7.4
Storage	Store at -20...-70°C upon receipt. Avoid multiple freeze-thaw cycles
Note	For research use only
Isotype	mIgG1
Clonality	Recombinant
Uniprot ID	NA
Expiration Date	12 months from date of receipt.