



Product Datasheet

Mouse monoclonal antibody to HPV-18 E7 (orb1108649)

Description	Mouse monoclonal antibody to HPV-18 E7
Species/Host	Mouse
Conjugation	Unconjugated
Tested Applications	ELISA, IF, IP, WB
Immunogen	Recombinant full-length HPV-18 E7-protein purified from E. coli
Preservatives	PBS pH 7.4, with 0.1% sodium azide
Storage	Store at -20...-70°C upon receipt. Divide antibody into aliquots prior usage. Avoid multiple freeze-thaw cycles.
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
Dilution Range	ELISA 0.05 – 0.1 µg/ml ; WB 1 µg/ml; IF 0.33 – 10 µg/ml; IP 2 µg/ml. Monoclonal antibody working titer has to be established practically for each particular antigen and assay format
Expiration Date	12 months from date of receipt.



IF analysis of HPV-18 E7 protein in U2OS



Western-Blot analysis of HPV-18 E7 prote...