

Product Datasheet

TLR6 antibody (orb1108873)

Description	Mouse monoclonal antibody to TLR6
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FC, Functional Studies, IF, IP, PS
Form/Appearance	1 ml (100 µg/ml) 0.2 µm filtered antibody solution in PBS, containing 0.1% bovine serum albumin.
Storage	Product should be stored at 4degC. Under recommended storage conditions, product is stable for one year.
Note	For research use only
Isotype	Mouse IgG1
Clonality	Monoclonal
Expiration Date	12 months from date of receipt.