

Product Datasheet

C5/C5a antibody (orb1108985)

Description	Mouse monoclonal antibody to C5/C5a
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	Functional Studies, IA, WB
Form/Appearance	1 ml (100 µg/ml) 0.2 µm filtered protein G purified antibody solution in PBS containing 0.1% bovine serum albumin. The endotoxin concentration is 24 EU/mg, determined with HIT302 LAL Assay.
Storage	Product should be stored at 4degC. Under recommended storage conditions, product is stable for one year.
Note	For research use only
Application notes	IA: antibody was used as detection antibody in ELISA. (Ref.1) W: used as ascites at a 1/1000 dilution. Incubation 2h on nitrocellulose blotted samples. (Ref.1)
Isotype	Mouse IgG2a
Clonality	Monoclonal
Expiration Date	12 months from date of receipt.