

## Product Datasheet

### Nectin-2 antibody (FITC) (orb1109147)

<b>Description</b>	Mouse monoclonal antibody against Nectin-2 conjugated to FITC
<b>Species/Host</b>	Mouse
<b>Reactivity</b>	Mouse
<b>Conjugation</b>	FITC
<b>Tested Applications</b>	IA, IF, IHC-Fr, IP, WB
<b>Form/Appearance</b>	1 ml (100 µg/ml) 0.2 µm protein G purified FITC conjugated antibody solution in PBS, containing 1% bovine serum albumin and 0.02% sodium azide.
<b>Storage</b>	Product should be stored at 4degC. Under recommended storage conditions, product is stable for at least one year.
<b>Note</b>	For research use only
<b>Application notes</b>	F: Fixation with PFA or acetone/methanol is recommended. IA: HM1052 can be used as a detection antibody. W: proteins were separated by SDS-PAGE and electrophoretically transferred onto polyvinylidene fluoride sheets, which were then incubated with primary antibodies. Binding was detected with HRP-conjugated secondary antibodies and Western blotting detection reagents (Ref.3).
<b>Isotype</b>	Rat IgG2a
<b>Clonality</b>	Monoclonal
<b>Expiration Date</b>	12 months from date of receipt.