

Product Datasheet

MPO antibody (FITC) (orb1109149)

Description	Mouse monoclonal antibody against MPO conjugated to FITC
Species/Host	Mouse
Reactivity	Mouse
Conjugation	FITC
Tested Applications	FC, IA, IF, IHC-Fr
Form/Appearance	0.5 ml (100 µg/ml) 0.2 um filtered protein G purified FITC conjugated antibody solution in PBS, containing 1.0% bovine serum albumin and 0.02% sodium azide.
Storage	Product should be stored at 4degC. Under recommended storage conditions, product is stable for one year.
Note	For research use only
Application notes	IA: the antibody can be used as a detection antibody. FC: Cells are fixed with 4% paraformaldehyde and permeabilized with 0.1% saponin. IA: HM1051 was used as a capture antibody. The experiment performed was a specificity experiment with different concentrations of HM1051.
Isotype	Mouse IgG1
Clonality	Monoclonal
Expiration Date	12 months from date of receipt.