

Product Datasheet

MBL-C antibody (orb1109170)

Description	Mouse monoclonal antibody to MBL-C
Species/Host	Mouse
Reactivity	Mouse
Conjugation	Unconjugated
Tested Applications	IA, WB
Form/Appearance	1 ml (100 µg/ml) 0.2 µm filtered antibody solution in PBS, containing 0.02% sodium azide and 0.1% bovine serum albumin.
Storage	Product should be stored at 4degC. Under recommended storage conditions, product is stable for one year.
Note	For research use only
Application notes	The monoclonal antibody 16A8 can be used for immuno assay and Western blotting (non-reduced). The antibody 16A8 is not useful for immunohistology.
Isotype	Rat IgG2a
Clonality	Monoclonal
Expiration Date	12 months from date of receipt.