

Product Datasheet

PITPNM2 antibody (orb1124794)

Description	Rabbit polyclonal antibody to PITPNM2
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, WB
Immunogen	raised against a recombinant human PITPNM2 protein 162-435aa (BC141813)
Preservatives	PBS, pH 7.4 with 0.02% Sodium Azide
Form/Appearance	Liquid
Concentration	100ug/100ul
Storage	This product is stable for several weeks at 4°C as an undiluted liquid. Dilute only prior to immediate use. For extended storage, aliquot contents and freeze at -20°C or below. Avoid cycles of freezing and thawing. Expiration date is one (1) year from date of receipt.
Note	For research use only
Application notes	ELISA titer: 1:20,000-1:80,000; Western blot analysis: 1:500-1:1,000; Immunohistochemistry: 1:100-1:500; Optimal dilutions/concentrations should be determined by the end use
Clonality	Polyclonal
Dilution Range	ELISA titer: 1:20,000-1:80,000; Western blot analysis: 1:500-1:1,000; Immunohistochemistry: 1:100-1:500; Optimal dilutions/concentrations should be determined by the end use
Expiration Date	12 months from date of receipt.