

Product Datasheet

Trx Tag antibody (orb1125307)

Description	Mouse monoclonal antibody to Trx Tag
Species/Host	Mouse
Reactivity	Other
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Recombinant protein corresponding to fragment (aa/1-111) Trx Tag .
Preservatives	citrate-tris-HCl buffer pH 7.0 with 0.02% Proclin 300.
Form/Appearance	Liquid
Concentration	100ug/100ul
Storage	This product is stable for several weeks at 4° C as an undiluted liquid. Dilute only prior to immediate use. For extended storage, aliquot contents and freeze at -20° C or below. Avoid cycles of freezing and thawing. Expiration date is one (1) year from date of receipt.
Note	For research use only
Application notes	WB (1:1000) to Trx-tagged fusion protein with Anti-Trx Tag
Isotype	IgG1
Clonality	Monoclonal
Dilution Range	WB (1:1000) to Trx-tagged fusion protein with Anti-Trx Tag
Expiration Date	12 months from date of receipt.