

## Product Datasheet

### GFP Tag antibody (orb1125308)

<b>Description</b>	Mouse monoclonal antibody to GFP Tag
<b>Species/Host</b>	Mouse
<b>Reactivity</b>	Other
<b>Conjugation</b>	Unconjugated
<b>Tested Applications</b>	ELISA, ICC
<b>Immunogen</b>	Recombinant protein corresponding to full length GFP
<b>Preservatives</b>	citrate-tris-HCl buffer pH 7.0 with 0.02% Proclin 300.
<b>Form/Appearance</b>	Liquid
<b>Concentration</b>	100ug/100ul
<b>Storage</b>	This product is stable for several weeks at 4° C as an undiluted liquid. Dilute only prior to immediate use. For extended storage, aliquot contents and freeze at -20° C or below. Avoid cycles of freezing and thawing. Expiration date is one (1) year from date of receipt.
<b>Note</b>	For research use only
<b>Application notes</b>	1. Image taken under influorescence microscope;2. ICC detection (1:1);3. DAPI staining; 4. Merge image 1-3.
<b>Isotype</b>	IgG1
<b>Clonality</b>	Monoclonal
<b>Dilution Range</b>	1. Image taken under influorescence microscope;2. ICC detection (1:1);3. DAPI staining; 4. Merge image 1-3.
<b>Expiration Date</b>	12 months from date of receipt.