

Product Datasheet

Ecdysone Receptor (EcR) B1 (orb1129156)

Description	Mouse monoclonal antibody to Ecdysone Receptor (EcR) B1
Species/Host	Mouse
Reactivity	Insect
Conjugation	Unconjugated
Preservatives	PBS, pH 7.4 with 0.02% Sodium Azide
Form/Appearance	Liquid
Storage	This product is stable for several weeks at 4° C as an undiluted liquid. Dilute only prior to immediate use. For extended storage, aliquot contents and freeze at -20° C or below. Avoid cycles of freezing and thawing. Expiration date is one (1) year from date of receipt.
Note	For research use only
Isotype	IgG2b
Clonality	Monoclonal
Expiration Date	12 months from date of receipt.