

Product Datasheet

Anti-ZMAT3 (orb1132938)

Description	Rabbit polyclonal antibody to ZMAT3
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Raised against a recombinate human ZMAT3 protein 1-289aa (BC002896).
Preservatives	PBS, pH 7.4 with 0.02% Sodium Azide
Form/Appearance	Liquid
Concentration	100ug/100ul
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Application notes	ELISA titer: 1:20,000-1:80,000 Western blot analysis: 1:500-1:1,000. Optimal dilutions/concentrations should be determined by the end use.
Clonality	Polyclonal
Purity	Affinity purified
Dilution Range	ELISA titer: 1:20,000-1:80,000 Western blot analysis: 1:500-1:1,000. Optimal dilutions/concentrations should be determined by the end use.
Expiration Date	12 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann's Road,
Cambridge, CB5 8LA, United Kingdom
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\) 1223 859-353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)

Biorbyt LLC.

68 TW Alexander Drive,
Durham, NC, 27713, United States
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)9065211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)