

Product Datasheet

Anti-ACTL8 (orb1133305)

Description	Rabbit polyclonal antibody to ACTL8
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC
Immunogen	Rabbit polyclonal ACTL8 antibody was raised against a peptide near intermediate residues of human ACTL8 protein.
Preservatives	PBS, pH 7.4 with 0.02% Sodium Azide
Form/Appearance	Liquid
Concentration	100ug/100ul
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Application notes	Optimal dilutions to be determined by researchers
Isotype	IgG
Clonality	Polyclonal
Antibody Type	Primary Antibody
Purity	Affinity purified
Dilution Range	Optimal dilutions to be determined by researchers
Expiration Date	12 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann's Road,
Cambridge, CB5 8LA, United Kingdom
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\) 1223 859-353](tel:+44201223859353) | Fax: [+1 \(415\) 651-8558](tel:+14156518558)

Biorbyt LLC.

68 TW Alexander Drive,
Durham, NC, 27713, United States
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+14159065211) | Fax: [+1 \(415\) 651-8558](tel:+14156518558)