

---

## Product Datasheet

### S-Tag-UNLB antibody (orb1133730)

**Description**

Mouse monoclonal antibody to S-Tag-UNLB

**Species/Host**

Mouse

**Reactivity**

Other

**Conjugation**

Unconjugated

**Tested Applications**

FC

**Immunogen**

S-Tag recombinant protein

**Preservatives**

100 mM borate buffered saline, pH 8.0. No preservatives or amine-containing buffer salts added.

**Form/Appearance**

Liquid

**Concentration**

100ug/200ul

**Storage**

This product is stable for several weeks at 4° C as an undiluted liquid. Dilute only prior to immediate use. For extended storage, aliquot contents and freeze at -20° C or below. Avoid cycles of freezing and thawing. Expiration date is one (1) year from date of receipt.

**Note**

For research use only

**Application notes**

• Flow cytometry • Enzyme linked immunosorbent assay (ELISA); Clone SBSTAGa can be paired with SBSTAGb or SBSTAGc. • Protein purification

**Isotype**

Mouse IgG1k

**Clonality**

Monoclonal

**Dilution Range**

• Flow cytometry • Enzyme linked immunosorbent assay (ELISA); Clone SBSTAGa can be paired with SBSTAGb or SBSTAGc. • Protein purification

**Expiration Date**

12 months from date of receipt.