


Product Datasheet

P105RB (phospho-Ser780) antibody (orb11338)

| | |
|----------------------------|---|
| Description | P105RB (phospho-Ser780) antibody |
| Species/Host | Rabbit |
| Reactivity | Human, Mouse, Rat |
| Conjugation | Unconjugated |
| Tested Applications | FC, ICC, IHC-P, WB |
| Immunogen | KLH conjugated Synthesised phosphopeptide derived from human Rb around the phosphorylation site of Ser780 (L(p-S)PI) |
| Target | phospho-Rb (Ser780) |
| Preservatives | 0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol. |
| Form/Appearance | Liquid |
| Concentration | 1mg/ml |
| Storage | Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles. |
| Note | For research use only |
| Isotype | IgG |
| Clonality | Polyclonal |
| MW | 106 kDa |
| Uniprot ID | P06400 |
| Dilution Range | WB=1:500-2000, IHC-P=1:100-500, ICC=1:100, Flow-Cyt=1ug/test, IHC-F=1:100-500, IF=1:100-500, ELISA=1:5000-10000 |
| Expiration Date | 12 months from date of receipt. |



IHC-P analysis of mouse colon using phos...