

Product Datasheet

Influenza A virus HA tag antibody (orb1134542)

Description	Rabbit polyclonal antibody to Influenza A virus HA tag
Species/Host	Rabbit
Reactivity	Other
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Peptide
Preservatives	PBS, pH 7.4 with 0.1% sodium azide, 50% Glycerol
Form/Appearance	Liquid
Concentration	1ug/ul
Storage	This product is stable for several weeks at 4° C as an undiluted liquid. Dilute only prior to immediate use. For extended storage, aliquot contents and freeze at -20° C or below. Avoid cycles of freezing and thawing. Expiration date is one (1) year from date of receipt.
Note	For research use only
Application notes	ELISA 1:1,000~100,000 Western blotting 1:100~1,000
Isotype	IgG
Clonality	Polyclonal
Dilution Range	ELISA 1:1,000~100,000 Western blotting 1:100~1,000
Expiration Date	12 months from date of receipt.