



Product Datasheet

Sex Hormone Binding Globulin antibody (orb11366)

Description	Sex Hormone Binding Globulin antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IF, IHC-Fr, IHC-P, WB
Immunogen	KLH conjugated synthetic peptide derived from human SHBG (301-402/402 aa)
Target	SHBG
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	44 kDa
Uniprot ID	A0A6P7EHR7
Expiration Date	12 months from date of receipt.



Western blot analysis of mouse testis ly...



Paraffin-embedded rat uterus tissue fixe...