
Product Datasheet

Syk (Ab-323) antibody (orb1139570)

Description

Rabbit polyclonal antibody to Syk (Ab-323)

Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC, IP, WB
Immunogen	The antibody was produced against synthesized non-phosphopeptide derived from human syk around the phosphorylation site of tyrosine 323 (N-P-YP-E-P).
Preservatives	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Form/Appearance	Liquid
Concentration	100ug/100ul
Storage	This product is stable for several weeks at 4° C as an undiluted liquid. Dilute only prior to immediate use. For extended storage, aliquot contents and freeze at -20° C or below. Avoid cycles of freezing and thawing. Expiration date is one (1) year from date of receipt.
Note	For research use only
Application notes	WB: 1:500~1:1000 IHC: 1:50~1:100 IF:1:100~1:200
Isotype	IgG
Clonality	Polyclonal
Dilution Range	WB: 1:500~1:1000 IHC: 1:50~1:100 IF:1:100~1:200
Expiration Date	12 months from date of receipt.