

Product Datasheet

Smad3 (Ab-425) antibody (orb1139605)

Description

Rabbit polyclonal antibody to Smad3 (Ab-

Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IHC, IP
Immunogen	The antibody was produced against synthesized non-phosphopeptide derived from human Smad3 around the phosphorylation site of serine 425 (C-S-S-V-SP).
Preservatives	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Form/Appearance	Liquid
Concentration	100ug/100ul
Storage	This product is stable for several weeks at 4° C as an undiluted liquid. Dilute only prior to immediate use. For extended storage, aliquot contents and freeze at -20° C or below. Avoid cycles of freezing and thawing. Expiration date is one (1) year from date of receipt.
Note	For research use only
Application notes	IHC: 1:50~1:100 IF:1:100~1:200
Isotype	IgG
Clonality	Polyclonal
Dilution Range	IHC: 1:50~1:100 IF:1:100~1:200
Expiration Date	12 months from date of receipt.