
Product Datasheet

IRS-1 (Ab-639) antibody (orb1139691)

Description

Rabbit polyclonal antibody to IRS-1 (Ab-639)

Species/Host

Rabbit

Reactivity

Human, Mouse, Rat

Conjugation

Unconjugated

Tested Applications

IP, WB

Immunogen

The antibody was produced against synthesized non-phosphopeptide derived from human IRS-1 around the phosphorylation site of serine 639 (P-K-SP-V-S).

Preservatives

Rabbit IgG in phosphate buffered saline (without Mg^{2+} and Ca^{2+}), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Form/Appearance

Liquid

Concentration

100ug/100ul

Storage

This product is stable for several weeks at 4° C as an undiluted liquid. Dilute only prior to immediate use. For extended storage, aliquot contents and freeze at -20° C or below. Avoid cycles of freezing and thawing. Expiration date is one (1) year from date of receipt.

Note

For research use only

Application notes

WB: 1:500~1:1000 IF:1:100~1:200

Isotype

IgG

Clonality

Polyclonal

Dilution Range

WB: 1:500~1:1000 IF:1:100~1:200

Expiration Date

12 months from date of receipt.