
Product Datasheet

Cyclin B1 (phospho-Ser147) antibody (orb1139912)

Description

Rabbit polyclonal antibody to Cyclin B1

Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ICC, IHC
Immunogen	The antibody was produced against synthesized phosphopeptide derived from human Cyclin B1 around the phosphorylation site of Serine 147 (A-F-SP-D-V).
Preservatives	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Form/Appearance	Liquid
Concentration	100ug/100ul
Storage	This product is stable for several weeks at 4° C as an undiluted liquid. Dilute only prior to immediate use. For extended storage, aliquot contents and freeze at -20° C or below. Avoid cycles of freezing and thawing. Expiration date is one (1) year from date of receipt.
Note	For research use only
Application notes	IHC: 1:50~100 IF:1:100~1:200
Isotype	IgG
Clonality	Polyclonal
Dilution Range	IHC: 1:50~100 IF:1:100~1:200
Expiration Date	12 months from date of receipt.