
Product Datasheet

PLC-gamma 2(phospho-Tyr1217) antibody (orb1139918)

Description

Rabbit polyclonal antibody to PLC-

Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IHC
Immunogen	The antibody was produced against synthesized phosphopeptide derived from human PLC- γ 2 around the phosphorylation site of tyrosine 1217 (F-L-YP-D-T).
Preservatives	PBS, pH 7.4 with 0.02% sodium azide.
Form/Appearance	Liquid
Concentration	100ug/100ul
Storage	This product is stable for several weeks at 4° C as an undiluted liquid. Dilute only prior to immediate use. For extended storage, aliquot contents and freeze at -20° C or below. Avoid cycles of freezing and thawing. Expiration date is one (1) year from date of receipt.
Note	For research use only
Application notes	IHC: 1:50~100
Isotype	IgG
Clonality	Polyclonal
Dilution Range	IHC: 1:50~100
Expiration Date	12 months from date of receipt.