
Product Datasheet

**HDAC4/HDAC5/HDAC9 (phospho-Ser246/259/220)
antibody (orb1139921)**

Description

Rabbit polyclonal antibody to

Ser246/259/220)

Species/Host

Rabbit

Reactivity

Human, Mouse

Conjugation

Unconjugated

Tested Applications

IHC, WB

Immunogen

The antibody was produced against synthesized phosphopeptide derived from human HDAC4/HDAC5/HDAC9 around the phosphorylation site of serine 246/259/220 (T-A-SP-EP).

Preservatives

Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Form/Appearance

Liquid

Concentration

100ug/100ul

Storage

This product is stable for several weeks at 4° C as an undiluted liquid. Dilute only prior to immediate use. For extended storage, aliquot contents and freeze at -20° C or below. Avoid cycles of freezing and thawing. Expiration date is one (1) year from date of receipt.

Note

For research use only

Application notes

WB: 1:500~1:1000 IHC:1:50~1:100

Isotype

IgG

Clonality

Polyclonal

Dilution Range

WB: 1:500~1:1000 IHC:1:50~1:100

Expiration Date

12 months from date of receipt.