
Product Datasheet

eIF4B (phospho-Ser422) antibody (orb1139923)

Description

Rabbit polyclonal antibody to eIF4B

Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IHC, WB
Immunogen	The antibody was produced against synthesized phosphopeptide derived from human eIF4B around the phosphorylation site of Serine 422 (T-G-SP-E-S).
Preservatives	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Form/Appearance	Liquid
Concentration	100ug/100ul
Storage	This product is stable for several weeks at 4° C as an undiluted liquid. Dilute only prior to immediate use. For extended storage, aliquot contents and freeze at -20° C or below. Avoid cycles of freezing and thawing. Expiration date is one (1) year from date of receipt.
Note	For research use only
Application notes	WB: 1:500~1000 IHC: 1:50~100
Isotype	IgG
Clonality	Polyclonal
Dilution Range	WB: 1:500~1000 IHC: 1:50~100
Expiration Date	12 months from date of receipt.