

Product Datasheet

SHP-2 (Phospho-Tyr580) antibody (orb1139942)

| | |
|----------------------------|--|
| Description | Rabbit polyclonal antibody to SHP-2 (Phospho-Tyr580) |
| Species/Host | Rabbit |
| Reactivity | Human, Mouse, Rat |
| Conjugation | Unconjugated |
| Tested Applications | IHC |
| Immunogen | The antiserum was produced against synthesized phosphopeptide derived from human SHP-2 around the phosphorylation site of tyrosine 580 (R-V-YP-E-N). |
| Preservatives | phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol. |
| Form/Appearance | Liquid |
| Concentration | 1ug/ul |
| Storage | This product is stable for several weeks at 4° C as an undiluted liquid. Dilute only prior to immediate use. For extended storage, aliquot contents and freeze at -20° C or below. Avoid cycles of freezing and thawing. Expiration date is one (1) year from date of receipt. |
| Note | For research use only |
| Application notes | IHC: 1:50-1:100 |
| Clonality | Polyclonal |
| Dilution Range | IHC: 1:50-1:100 |
| Expiration Date | 12 months from date of receipt. |