
Product Datasheet

**NF-Kappa B p65 (phospho-Ser311) antibody
(orb1139987)**

Description

Rabbit polyclonal antibody to NF-

Species/Host

Rabbit

Reactivity

Human, Mouse, Rat

Conjugation

Unconjugated

Tested Applications

IHC

Immunogen

The antiserum was produced against synthesized phosphopeptide derived from human NFκB-p65 around the phosphorylation site of serine 311 (F-K-SP-I-M).

Preservatives

phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Form/Appearance

Liquid

Concentration

1ug/ul

Storage

This product is stable for several weeks at 4° C as an undiluted liquid. Dilute only prior to immediate use. For extended storage, aliquot contents and freeze at -20° C or below. Avoid cycles of freezing and thawing. Expiration date is one (1) year from date of receipt.

Note

For research use only

Application notes

IHC: 1:50~1:100

Clonality

Polyclonal

Dilution Range

IHC: 1:50~1:100

Expiration Date

12 months from date of receipt.