
Product Datasheet

SAPK/JNK (Phospho-Thr183) antibody (orb1139996)

Description

Rabbit polyclonal antibody to SAPK/JNK

Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IHC
Immunogen	The antiserum was produced against synthesized phosphopeptide derived from human JNK2/3/SAPK around the phosphorylation site of threonine 183 (M-M-TP-P-Y).
Preservatives	phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Form/Appearance	Liquid
Concentration	1ug/ul
Storage	This product is stable for several weeks at 4° C as an undiluted liquid. Dilute only prior to immediate use. For extended storage, aliquot contents and freeze at -20° C or below. Avoid cycles of freezing and thawing. Expiration date is one (1) year from date of receipt.
Note	For research use only
Application notes	IHC: 1:50~1:100
Clonality	Polyclonal
Dilution Range	IHC: 1:50~1:100
Expiration Date	12 months from date of receipt.