

Product Datasheet

c-Kit (Phospho-Tyr721) antibody (orb1140003)

Description	Rabbit polyclonal antibody to c-Kit (Phospho-Tyr721)
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	The antiserum was produced against synthesized phosphopeptide derived from human c-kit around the phosphorylation site of tyrosine 721 (N-E-YP-M-D).
Preservatives	phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Form/Appearance	Liquid
Concentration	1ug/ul
Storage	This product is stable for several weeks at 4° C as an undiluted liquid. Dilute only prior to immediate use. For extended storage, aliquot contents and freeze at -20° C or below. Avoid cycles of freezing and thawing. Expiration date is one (1) year from date of receipt.
Note	For research use only
Application notes	WB: 1:500~1:1000
Clonality	Polyclonal
Dilution Range	WB: 1:500~1:1000
Expiration Date	12 months from date of receipt.