

---

## Product Datasheet

### Tyrosine Hydroxylase (Phospho-Ser40) antibody (orb1140027)

**Description**

Rabbit polyclonal antibody to

Ser40)

**Species/Host**

Rabbit

**Reactivity**

Human, Mouse, Rat

**Conjugation**

Unconjugated

**Tested Applications**

ICC

**Immunogen**

The antibody was produced against synthesized phosphopeptide derived from human Tyrosine Hydroxylase around the phosphorylation site of serine 40(R-Q-SP-L-I).

**Preservatives**

Rabbit IgG in phosphate buffered saline (without Mg<sup>2+</sup> and Ca<sup>2+</sup>), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

**Form/Appearance**

Liquid

**Concentration**

100ug/100ul

**Storage**

This product is stable for several weeks at 4° C as an undiluted liquid. Dilute only prior to immediate use. For extended storage, aliquot contents and freeze at -20° C or below. Avoid cycles of freezing and thawing. Expiration date is one (1) year from date of receipt.

**Note**

For research use only

**Application notes**

IF:1:100~1:200

**Isotype**

IgG

**Clonality**

Polyclonal

**Dilution Range**

IF:1:100~1:200

**Expiration Date**

12 months from date of receipt.