
Product Datasheet

**IκB-alpha (Phospho-Ser32/Ser36) antibody
(orb1140075)**

Description

Rabbit polyclonal antibody to IκB-

Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IHC
Immunogen	The antiserum was produced against synthesized phosphopeptide derived from human IκB-α around the phosphorylation site of serine 32/36 (H-D-SP-G-L-D-SP -M-K).
Preservatives	phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Form/Appearance	Liquid
Concentration	1ug/ul
Storage	This product is stable for several weeks at 4° C as an undiluted liquid. Dilute only prior to immediate use. For extended storage, aliquot contents and freeze at -20° C or below. Avoid cycles of freezing and thawing. Expiration date is one (1) year from date of receipt.
Note	For research use only
Application notes	IHC: 1:50~1:100
Clonality	Polyclonal
Dilution Range	IHC: 1:50~1:100
Expiration Date	12 months from date of receipt.