

---

## Product Datasheet

**beta-Catenin (Phospho-Thr41/Ser45) antibody  
(orb1140108)**

**Description**

Rabbit polyclonal antibody to beta-

**Species/Host**

Rabbit

**Reactivity**

Human, Mouse, Rat

**Conjugation**

Unconjugated

**Tested Applications**

WB

**Immunogen**

The antiserum was produced against synthesized phosphopeptide derived from human $\beta$ -Catenin around the phosphorylation site of threonine 41/serine 45 (A-T-TP-T-A-P-SP-L-S).

**Preservatives**

phosphate buffered saline (without Mg<sup>2+</sup> and Ca<sup>2+</sup>), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

**Form/Appearance**

Liquid

**Concentration**1 $\mu$ g/ $\mu$ l**Storage**

This product is stable for several weeks at 4° C as an undiluted liquid. Dilute only prior to immediate use. For extended storage, aliquot contents and freeze at -20° C or below. Avoid cycles of freezing and thawing. Expiration date is one (1) year from date of receipt.

**Note**

For research use only

**Application notes**

WB: 1:500~1:1000

**Clonality**

Polyclonal

**Dilution Range**

WB: 1:500~1:1000

**Expiration Date**

12 months from date of receipt.