
Product Datasheet

Tau (Phospho-Ser356) antibody (orb1140123)

Description

Rabbit polyclonal antibody to Tau

Species/Host

Rabbit

Reactivity

Human, Mouse, Rat

Conjugation

Unconjugated

Tested Applications

ICC, WB

Immunogen

The antibody was produced against synthesized phosphopeptide derived from human Tau around the phosphorylation site of serine 356 (I-G-SP-L-D).

Preservatives

Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Form/Appearance

Liquid

Concentration

100ug/100ul

Storage

This product is stable for several weeks at 4° C as an undiluted liquid. Dilute only prior to immediate use. For extended storage, aliquot contents and freeze at -20° C or below. Avoid cycles of freezing and thawing. Expiration date is one (1) year from date of receipt.

Note

For research use only

Application notes

WB: 1:500~1:1000 IF:1:100~1:200

Isotype

IgG

Clonality

Polyclonal

Dilution Range

WB: 1:500~1:1000 IF:1:100~1:200

Expiration Date

12 months from date of receipt.