

Product Datasheet

V5 Tag antibody (orb1140356)

Description	Rabbit polyclonal antibody to V5 Tag
Species/Host	Rabbit
Reactivity	Other
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	GKPIPPLLGLDST peptide.
Form/Appearance	Liquid
Concentration	100ug/100ul
Storage	This product is stable for several weeks at 4° C as an undiluted liquid. Dilute only prior to immediate use. For extended storage, aliquot contents and freeze at -20° C or below. Avoid cycles of freezing and thawing. Expiration date is one (1) year from date of receipt.
Note	For research use only
Application notes	Starting dilution of 1:1000 for Western Blot and 1:10,000 for ELISA. Other applications may be possible and will need to be evaluated and optimized by end user.
Clonality	Polyclonal
Dilution Range	Starting dilution of 1:1000 for Western Blot and 1:10,000 for ELISA. Other applications may be possible and will need to be evaluated and optimized by end user.
Expiration Date	12 months from date of receipt.