

---

## Product Datasheet

### D Tag antibody (orb1140363)

**Description**

Mouse monoclonal antibody to D Tag

**Species/Host**

Mouse

**Reactivity**

Other

**Conjugation**

Unconjugated

**Tested Applications**

IP, WB

**Immunogen**

A synthetic peptide (DYKDDDDK) coupled to KLH.

**Form/Appearance**

Liquid

**Concentration**

100µg / 200µl

**Storage**

This product is stable for several weeks at 4° C as an undiluted liquid. Dilute only prior to immediate use. For extended storage, aliquot contents and freeze at -20° C or below. Avoid cycles of freezing and thawing. Expiration date is one (1) year from date of receipt.

**Note**

For research use only

**Application notes**

Optimal working dilutions should be determined experimentally by the investigator. Prepare working dilution immediately before use. Suggested starting dilutions are as follows: Western Blot at dilutions of 1:1000. Immunoprecipitation and immunostaining at dilutions of 1:1000.

**Isotype**

IgG2b with lambda light chain

**Clonality**

Monoclonal

**Dilution Range**

Optimal working dilutions should be determined experimentally by the investigator. Prepare working dilution immediately before use. Suggested starting dilutions are as follows: Western Blot at dilutions of 1:1000. Immunoprecipitation and immunostaining at dilutions of 1:1000.

**Expiration Date**

12 months from date of receipt.