

Product Datasheet

V5 Tag antibody (orb1140365)

Description	Mouse monoclonal antibody to V5 Tag
Species/Host	Mouse
Reactivity	Other
Conjugation	Unconjugated
Tested Applications	ELISA, ICC, IHC, WB
Immunogen	Anti-V5 is a mouse monoclonal antibody, affinity purified from culture supernatant of hybridoma cells grown in a bioreactor.
Form/Appearance	Liquid
Concentration	100ug/100ul
Storage	This product is stable for several weeks at 4° C as an undiluted liquid. Dilute only prior to immediate use. For extended storage, aliquot contents and freeze at -20° C or below. Avoid cycles of freezing and thawing. Expiration date is one (1) year from date of receipt.
Note	For research use only
Application notes	Western blot analysis: 1:1000 ELISA: 1:5000-10,000 IHC: 1:50-1:100. Other applications will need optimization from end users.
Isotype	IgG1
Clonality	Monoclonal
Dilution Range	Western blot analysis: 1:1000 ELISA: 1:5000-10,000 IHC: 1:50-1:100. Other applications will need optimization from end users.
Expiration Date	12 months from date of receipt.