

Product Datasheet

RFP Tag antibody (orb1140370)

Description	Rabbit polyclonal antibody to RFP Tag
Species/Host	Rabbit
Reactivity	Other
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, WB
Immunogen	The immunogen is a Red Fluorescent Protein (RFP) fusion protein corresponding to the full length amino acid sequence (234aa) derived from the mushroom polyp coral <i>Discosoma</i> .
Form/Appearance	Liquid
Concentration	100ug/100ul
Storage	This product is stable for several weeks at 4° C as an undiluted liquid. Dilute only prior to immediate use. For extended storage, aliquot contents and freeze at -20° C or below. Avoid cycles of freezing and thawing. Expiration date is one (1) year from date of receipt.
Note	For research use only
Application notes	ELISA 1:200,000; WESTERN BLOT 1:1,000; IMMUNOHISTOCHEMISTRY 1:1,000
Isotype	IgG
Clonality	Polyclonal
Dilution Range	ELISA 1:200,000; WESTERN BLOT 1:1,000; IMMUNOHISTOCHEMISTRY 1:1,000
Expiration Date	12 months from date of receipt.