

Product Datasheet

HA Tag antibody (orb1140372)

Description	Rabbit polyclonal antibody to HA Tag
Species/Host	Rabbit
Reactivity	Other
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	The antiserum was produced by immunizing rabbit with synthesized peptide containing the influenza hemagglutinin epitope (Y-P-Y-D-V-P-D-Y-A).
Form/Appearance	Liquid
Concentration	100ug/100ul
Storage	This product is stable for several weeks at 4° C as an undiluted liquid. Dilute only prior to immediate use. For extended storage, aliquot contents and freeze at -20° C or below. Avoid cycles of freezing and thawing. Expiration date is one (1) year from date of receipt.
Note	For research use only
Application notes	WB: 1:500~1:1000
Isotype	IgG
Clonality	Polyclonal
Dilution Range	WB: 1:500~1:1000
Expiration Date	12 months from date of receipt.