

Product Datasheet

His Tag antibody (orb1140375)

Description	Rabbit polyclonal antibody to His Tag
Species/Host	Rabbit
Reactivity	Other
Conjugation	Unconjugated
Tested Applications	ELISA, ICC, IHC, IP, WB
Immunogen	The antiserum was produced by immunizing rabbit with synthesized peptide containing the His epitope (H-H-H-H-H).
Form/Appearance	Liquid
Concentration	100ug/100ul
Storage	This product is stable for several weeks at 4°C as an undiluted liquid. Dilute only prior to immediate use. For extended storage, aliquot contents and freeze at -20°C or below. Avoid cycles of freezing and thawing. Expiration date is one (1) year from date of receipt.
Note	For research use only
Application notes	ELISA 1: 2500; WB 1:500-1000; IP 5 µg/sample; IF 5µg/mL; IHC 5 µg/mL
Isotype	IgG
Clonality	Polyclonal
Dilution Range	ELISA 1: 2500; WB 1:500-1000; IP 5 µg/sample; IF 5µg/mL; IHC 5 µg/mL
Expiration Date	12 months from date of receipt.